


# EVERYTHING YOU NEED TO SWIFTLY TURNOVER THE OPERATING ROOM. SIMPLIFIED.



# Improve Efficiencies in the Operating Room with Quick, Safe, Cleaning Solutions. Simplified by HALYARD\*

Add some simplicity and efficiency back into your workday with the complete line of HALYARD\* Room Turnover Kits. Featuring products you've known and trusted for years, our standard and custom kits make it easier to get everything you need for a clean, safe operating room (OR) in your acute or outpatient surgical facility.

This family of high-quality products, brought to you by the makers of AERO SERIES\* Surgical Gowns, delivers convenience at a cost-effective price.

Before, during and after surgery, HALYARD\* products support patient and staff safety in the OR.

---

## EXPLORE OUR POPULAR TURNOVER KIT PRODUCTS

Available for purchase in standard kits, this assortment of HALYARD\* disposable linens, helps staff members prevent infections in the OR and protects equipment from damage.



68474

### TABLE COVERS

---

- Quilted, quick absorbing table covers safeguard patients and equipment from fluids.
- The breathable, top layer of material provides a comfortable surface for direct patient contact.
- Available in two different finishes, covers trap fluids in the inner core for ease of clean up following procedures.
- The low linting design minimizes wrinkling to ensure a dry outer layer for protecting both patient skin and surrounding equipment.
- The impervious bottom layer helps prevent strikethrough and cross-contamination.
- Single use linens reduce the potential for biological and environmental contaminants being introduced in the OR.



### LIFT SHEET

---

Built for strength, the lift sheet supports up to 270 kg of weight.



### ARMBOARD COVERS

---

The fabric armboard covers with straps stay in place for added patient safety.



### HEADREST COVERS

---

With an impervious backing, headrest covers help prevent bacteria from colonizing.

## STANDARD TURNOVER KIT OPTIONS

### *Save Time Searching for Single-Use Essentials*

Easily pull a standard kit for each case and add it to your case cart for efficient operating room prep. Choose from multiple, space-saving kit options in convenient packs.


#### Room Turnover Kit I

| Code  | Components   | Quantity | Units Per Case |
|-------|--|----------|----------------|
| 68500 | Absorbent Table Cover, 101.6 cm x 228.6 cm, Grippy | 1        | 40             |
|       | Lift Sheet, 101.6 cm x 152.4 cm                    | 1        |                |




#### Room Turnover Kit II

| Code  | Components   | Quantity | Units Per Case |
|-------|--|----------|----------------|
| 68501 | Absorbent Table Cover, 101.6 cm x 228.6 cm, Grippy | 1        | 36             |
|       | Lift Sheet, 101.6 cm x 152.4 cm                    | 1        |                |
|       | Armboard Cover, 33 cm x 76.2 cm with Strap         | 2        |                |



#### Room Turnover Kit III

| Code  | Components   | Quantity | Units Per Case |
|-------|--|----------|----------------|
| 68502 | Absorbent Table Cover, 101.6 cm x 228.6 cm, Grippy | 1        | 36             |
|       | Lift Sheet, 101.6 cm x 152.4 cm                    | 1        |                |
|       | Armboard Cover, 33 cm x 76.2 cm with Strap         | 2        |                |
|       | Head Rest Cover, 35.5 cm x 38.1 cm                 | 1        |                |



## ENHANCE YOUR OR EFFICIENCY WITH OUR GRIPPY ABSORBENT TABLE COVER.

*This product can be ordered separately or as part of the following  
HALYARD\* Room Turnover Kits:*

- Room Turnover Kit I
- Room Turnover Kit II
- Room Turnover Kit III

| Code  | Description   | Units Per Case |
|-------|---|----------------|
| 68474 | Absorbent Table Cover, 101.6 cm x 228.6 cm, Grippy Finish | 20             |

Look to HALYARD\* for every essential product you need  
for smooth OR procedures and optimize your savings.  
From single-use products to capital equipment, HALYARD\*  
offers the right resources for safe and efficient patient care.





Jakelija:

# VINGMED

## Vingmed Oy

Olarinluoma 12 B, 02200 Espoo, Suomi  
+358 20 730 0350 | asiakaspalvelu@vingmed.fi

vingmed.fi



For more information,  
please contact your Halyard Customer Care Team:

UK & IE (toll-free)

Tel: 0800 085 82 04 - Fax: 0800 066 91 06  
[customerservice.uk.ie@owens-minor.com](mailto:customerservice.uk.ie@owens-minor.com)

Other countries

Tel: +32 2 585 51 14 - Fax: +32 2 897 00 95  
[customerservice.export@owens-minor.com](mailto:customerservice.export@owens-minor.com)

[www.halyardhealth.eu](http://www.halyardhealth.eu)



\*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2025. All rights reserved.  
008860 COPY-05307 - HC572-00 - EN GB